


euca

V_{CF} Line

Fast Exchange Semicompliant
PTCA Balloon Dilatation Catheter



eucatech

INTEGRATED MEDICINE

Balloon-Ø [mm]	1.5	2.0	2.5	2.75	3.0	3.25	3.5	4.0	4.5	5.0	5.5	6.0	6.5
Deflated [inch] Balloon Profile	0.022	0.023	0.024	0.024	0.025	0.025	0.026	0.027	0.029	0.030	0.032	0.034	0.036
Pressure [MPa/bar]	Compliance data [mm]												
0.6 / 6	1.46	1.92	2.42	2.64	2.84	3.12	3.38	3.80	4.33	4.86	5.30	5.70	6.30
0.7 / 7	1.48	1.96	2.46	2.69	2.92	3.18	3.44	3.90	4.41	4.92	5.40	5.85	6.40
0.8 / 8	1.50	2.00	2.50	2.75	3.00	3.25	3.50	4.00	4.50	5.00	5.50	6.00	6.50
0.9 / 9	1.53	2.02	2.52	2.79	3.06	3.30	3.54	4.10	4.55	5.05	5.60	6.10	6.60
1.0 / 10	1.55	2.04	2.54	2.83	3.13	3.35	3.58	4.20	4.60	5.10	5.65	6.18	6.70
1.1 / 11	1.57	2.06	2.56	2.87	3.19	3.40	3.62	4.30	4.65	5.15	5.70	6.26	6.80
1.2 / 12	1.59	2.08	2.58	2.90	3.22	3.44	3.66	4.35	4.70	5.20	5.75	6.32	6.90
1.3 / 13	1.61	2.09	2.60	2.92	3.25	3.47	3.70	4.39	4.73	5.23	5.80	6.37	6.97
1.4 / 14	1.63	2.10	2.62	2.95	3.28	3.50	3.74	4.43	4.76	5.26	5.85	6.42	7.04
1.5 / 15	1.65	2.11	2.64	2.97	3.30	3.53	3.78	4.47	4.79	5.29	5.90	6.47	7.11
1.6 / 16	1.67	2.12	2.66	3.00	3.33	3.56	3.82	4.50	4.82	5.32			
1.7 / 17	1.69	2.13	2.68	3.02	3.35	3.59	3.86	4.53					
1.8 / 18	1.71	2.14	2.70	3.04	3.37	3.62							
1.9 / 19	1.73	2.15	2.72	3.06									
2.0 / 20	1.75	2.16											

Nominal Pressure at 8 bar (0.8 MPa) ■
 Rated Burst Pressure ■



Order Information

Balloon Diameter	Order Number					
	Length					
	10 mm	15 mm	20 mm	25 mm	30 mm	40 mm
1.5 mm	euca V 1510	euca V 1515	euca V 1520	euca V 1525	-----	-----
2.0 mm	euca V 2010	euca V 2015	euca V 2020	euca V 2025	-----	-----
2.5 mm	euca V 2510	euca V 2515	euca V 2520	euca V 2525	euca V 2530	euca V 2540
2.75 mm	euca V 27510	euca V 27515	euca V 27520	euca V 27525	euca V 27530	euca V 27540
3.0 mm	euca V 3010	euca V 3015	euca V 3020	euca V 3025	euca V 3030	euca V 3040
3.25 mm	euca V 32510	euca V 32515	euca V 32520	euca V 32525	euca V 32530	euca V 32540
3.5 mm	euca V 3510	euca V 3515	euca V 3520	euca V 3525	euca V 3530	euca V 3540
4.0 mm	euca V 4010	euca V 4015	euca V 4020	euca V 4025	euca V 4030	euca V 4040
4.5 mm	euca V 4510	euca V 4515	euca V 4520	euca V 4525	euca V 4530	euca V 4540
5.0 mm	euca V 5010	euca V 5015	euca V 5020	euca V 5025	euca V 5030	euca V 5040
5.5 mm	euca V 5510	euca V 5515	euca V 5520	euca V 5525	euca V 5530	euca V 5540
6.0 mm	euca V 6010	euca V 6015	euca V 6020	euca V 6025	euca V 6030	euca V 6040
6.5 mm	euca V 6510	euca V 6515	euca V 6520	euca V 6525	euca V 6530	euca V 6540

Special balloon sizes on request, please contact us.



INTEGRATED MEDICINE

eucatech AG

Sieboldstr. 7 · D-97688 Bad Kissingen, Germany

Tel.: +49(0)971 7236-200 · Fax: +49(0)971 7236-201

www.eucatech.com · e-mail: eucatech@eucatech.de