



Camouflage®

Camouflage® Coating to mask the stent

Camouflage®
Coronary Stentsystem
with Athrombogenic Coating

 eucatech
INTEGRATED MEDICINE

Camouflage® Coronary Stentsystem with Athrombogenic Coating

Compliance Chart – Crimped Profile Camouflage®

Balloon-Ø [mm]	2.0	2.25	2.50	2.75	3.00	3.25	3.50	4.00	4.50	5.00	5.50	6.00
Profile [inch]	0.032	0.033	0.036	0.036	0.037	0.037	0.038	0.039	0.040	0.040	0.041	0.042
Pressure [MPa/bar]	Compliance data [mm]											
0.6 / 6	1.92	2.20	2.42	2.68	2.92	3.16	3.40	3.80	4.33	4.86	5.30	5.70
0.7 / 7	1.96	2.23	2.47	2.72	2.96	3.20	3.44	3.90	4.41	4.92	5.40	5.85
0.8 / 8	2.00	2.25	2.50	2.75	3.00	3.25	3.50	4.00	4.50	5.00	5.50	6.00
0.9 / 9	2.02	2.27	2.52	2.97	3.05	3.20	3.35	4.06	4.55	5.05	5.60	6.10
1.0 / 10	2.04	2.29	2.57	2.84	3.11	3.35	3.60	4.12	4.60	5.10	5.65	6.18
1.1 / 11	2.06	2.31	2.60	2.88	3.16	3.40	3.65	4.18	4.65	5.15	5.70	6.26
1.2 / 12	2.08	2.34	2.63	2.92	3.20	3.44	3.69	4.24	4.70	5.20	5.75	6.32
1.3 / 13	2.09	2.36	2.66	2.95	3.24	3.48	3.72	4.20	4.72	5.22	5.80	6.37
1.4 / 14	2.10	2.39	2.68	2.98	3.27	3.52	3.77	4.36	4.76	5.26	5.85	6.42
1.5 / 15	2.11	2.41	2.71	3.01	3.31	3.55	3.80	4.43	4.79	5.29	5.90	6.47
1.6 / 16	2.12	2.43	2.73	3.03	3.33	3.58	3.82	4.49	4.82	5.32		
1.7 / 17	2.13	2.44	2.74	3.05	3.35	3.62	3.87	4.55				
1.8 / 18	2.14	2.46	2.76	3.08	3.40	3.65	3.90	4.61				
1.9 / 19	2.15	2.48	2.79	3.11	3.44	3.69						
2.0 / 20	2.16	2.50	2.81									

■ Nominal Pressure 2.8 MPa ■ Nominal Pressure 3.5 MPa

Order Information Camouflage®

Small Vessels

Stent ID	Stent Length						
	8 mm	13 mm	18 mm	23 mm	28 mm	33 mm	38 mm
2.00 mm	CFL 20 / 8	CFL 20 / 13	CFL 20 / 18				
2.25 mm	CFL 225 / 8	CFL 225 / 13	CFL 225 / 18	CFL 225 / 23			
2.50 mm	CFL 25 / 8	CFL 25 / 13	CFL 25 / 18	CFL 25 / 23			

Standard Vessels

Stent ID	Stent Length						
	8 mm	13 mm	18 mm	23 mm	28 mm	33 mm	38 mm
2.75 mm	CFL 275 / 8	CFL 275 / 13	CFL 275 / 18	CFL 275 / 23	CFL 275 / 28		
3.00 mm	CFL 30 / 8	CFL 30 / 13	CFL 30 / 18	CFL 30 / 23	CFL 30 / 28	CFL 30 / 33	CFL 30 / 38
3.25 mm	CFL 325 / 8	CFL 325 / 13	CFL 325 / 18	CFL 325 / 23	CFL 325 / 28	CFL 325 / 33	CFL 325 / 38
3.50 mm	CFL 35 / 8	CFL 35 / 13	CFL 35 / 18	CFL 35 / 23	CFL 35 / 28	CFL 35 / 33	CFL 35 / 38
4.00 mm	CFL 40 / 8	CFL 40 / 13	CFL 40 / 18	CFL 40 / 23	CFL 40 / 28	CFL 40 / 33	CFL 40 / 38

Large Vessels

Stent ID	Stent Length						
	11 mm	13 mm	18 mm	23 mm	28 mm	33 mm	38 mm
4.50 mm	CFL 45 / 11	CFL 45 / 13	CFL 45 / 18	CFL 45 / 23	CFL 45 / 28	CFL 45 / 33	CFL 45 / 38
5.00 mm	CFL 50 / 11	CFL 50 / 13	CFL 50 / 18	CFL 50 / 23	CFL 50 / 28	CFL 50 / 33	CFL 50 / 38
5.50 mm	CFL 55 / 11	CFL 55 / 13	CFL 55 / 18	CFL 55 / 23	CFL 55 / 28	CFL 55 / 33	CFL 55 / 38
6.00 mm	CFL 60 / 11	CFL 60 / 13	CFL 60 / 18	CFL 60 / 23	CFL 60 / 28	CFL 60 / 33	CFL 60 / 38

Other sizes on request. Under our continuous R&D program, eucatech reserves the right to change product specifications without prior notice. Indications, contraindications, warnings and instructions for use can be found in the product box supplied with each device.

Label by Hoesly®



eucatech AG
Sieboldstrasse 7
D-97688 Bad Kissingen/Germany
Tel.: ++49(0)971-7236-200
Fax: ++49(0)971-7236-201
www.eucatech.com
E-Mail: eucatech@eucatech.de



INTEGRATED MEDICINE