


**euca**

# ***ST / STS* Line**

Coronary Stent - Stent System



 eucatech

INTEGRATED MEDICINE

Balloon-Ø [mm]	2.5	2.75	3.0	3.25	3.5	4.0
<b>Pressure [MPa/bar]</b>	<b>Compliance Data Stent Inner Diameter [mm]</b>					
0.2 / 2	2.42	2.66	2.92	3.17	3.42	3.92
0.3 / 3	2.44	2.68	2.94	3.18	3.44	3.94
0.4 / 4	2.46	2.71	2.96	3.21	3.46	3.96
0.5 / 5	2.48	2.73	2.98	3.23	3.48	3.98
<b>0.6 / 6</b>	<b>2.50</b>	<b>2.75</b>	<b>3.00</b>	<b>3.25</b>	<b>3.50</b>	<b>4.00</b>
0.7 / 7	2.51	2.77	3.02	3.26	3.52	4.03
0.8 / 8	2.52	2.79	3.06	3.27	3.56	4.04
0.9 / 9	2.54	2.81	3.08	3.29	3.59	4.09
1.0 / 10	2.56	2.84	3.11	3.31	3.62	4.11
1.1 / 11	2.59	2.86	3.13	3.33	3.65	4.13
1.2 / 12	2.61	2.88	3.16	3.35	3.68	4.17
1.3 / 13	2.64	2.90	3.19	3.37	3.71	4.22
1.4 / 14	2.66	2.92	3.22	3.40	3.74	4.26
1.5 / 15	2.69	2.94	3.26	3.43	3.77	4.29
1.6 / 16	2.72	2.96	3.29	3.45	3.80	4.34
1.7 / 17	2.75	2.99	3.33	3.47	3.83	4.39
1.8 / 18	2.79	3.00	3.36	3.49		
1.9 / 19	2.81	3.01	3.38			
2.0 / 20	2.83					

Nominal Pressure at (0.6 MPa) 6 bar; Rated Burst Pressure  



#### Order Information - Coronary Stent


Order Number							
Stent Length	8 mm	12 mm	16 mm	20 mm	24 mm	28 mm	32 mm
	euca ST 8	euca ST 12	euca ST 16	euca ST 20	euca ST 24	euca ST 28	euca ST 32

#### Order Information - Coronary Stent System

Order Number							
Stent I. D.	Stent Length						
	8 mm	12 mm	16 mm	20 mm	24 mm	28 mm	32 mm
2.50 mm	euca STS 2.5 / 8	euca STS 2.5 / 12	euca STS 2.5 / 16	euca STS 2.5 / 20			
2.75 mm	euca STS 2.75 / 8	euca STS 2.75 / 12	euca STS 2.75 / 16	euca STS 2.75 / 20	euca STS 2.75 / 24		
3.00 mm	euca STS 3.0 / 8	euca STS 3.0 / 12	euca STS 3.0 / 16	euca STS 3.0 / 20	euca STS 3.0 / 24	euca STS 3.0 / 28	euca STS 3.0 / 32
3.25 mm	euca STS 3.25 / 8	euca STS 3.25 / 12	euca STS 3.25 / 16	euca STS 3.25 / 20	euca STS 3.25 / 24	euca STS 3.25 / 28	euca STS 3.25 / 32
3.50 mm	euca STS 3.5 / 8	euca STS 3.5 / 12	euca STS 3.5 / 16	euca STS 3.5 / 20	euca STS 3.5 / 24	euca STS 3.5 / 28	euca STS 3.5 / 32
4.00 mm	euca STS 4.0 / 8	euca STS 4.0 / 12	euca STS 4.0 / 16	euca STS 4.0 / 20	euca STS 4.0 / 24	euca STS 4.0 / 28	euca STS 4.0 / 32



INTEGRATED MEDICINE



eucatech AG

Sieboldstr. 7 · D-97688 Bad Kissingen, Germany

Tel.: +49(0)971 7236-200 · Fax: +49(0)971 7236-201

www.eucatech.com · e-mail: eucatech@eucatech.de